

Our new BluPrint™ Dynamic Mixing Cartridge is really causing a stir.



Benco's best impression material is now available for your Dynamic Mixing Machine.

For a perfect void-free mix every time, rely on BluPrint Vinyl Polysiloxane (VPS). You'll save time and reduce remakes. Plus, you'll reduce the potential for cross-contamination. Our premium VPS impression material is thixotropic and fast-setting, with a three-dimensional resin and hydrophilic surfactant. It's compatible with gypsums, electroplating and refractory models. Sets in 1½ minutes. Vanilla mint flavor.

DYNAMIC MIXING MACHINE CARTRIDGE KIT

Contains: one 390ml cartridge and 10 mix tips.
 Heavy Body [3503-004]
 Monophase [3502-998]

HIGH PERFORMANCE CARTRIDGE KIT

Contains: two 50ml cartridges and 6 large, teal mix tips.
 Heavy Body [3179-544] Monophase [3179-562]
 Light Body [3179-571] Putty [3179-553]

ACCESSORY ITEMS

Dispensing Gun [2583-304]
 48 Large Teal Mix Tips [2530-352]
 100 Small Yellow Intraoral Tips [2530-361]

HANDMIXED PUTTY KIT [3179-535]

Contains: one base and catalyst (400g each) and two scoops.

Try These Other Fine BluPrint™ Products.

BLUPRINT BITE REGISTRATION MATERIAL

Fast, accurate impressions for on-the-spot repair models and lightning-quick provisional crowns. Our high-strength formula locks in the detail every time with no mess, no fuss. Available in regular, fast and super fast set.

HIGH PERFORMANCE CARTRIDGE KIT

Contains: two 76g (50ml) split cartridges and 6 mix tips.
 Regular Set [3180-121] - Sets in 2 minutes
 Fast Set [3180-130] - Sets in 1 minute
 Super Fast Set [3270-079] - Sets in 30 seconds

BLUPRINT BITE TRAYS

Simultaneously captures the master impression, counter-impression and bite registration. Ultra-fine MicroMesh™ webbing extends as the patient closes, allowing close inter-digitation of the teeth without contact.

BOX OF 35

Anterior [1746-827]
 Quadrant [1746-845]

BOX OF 50

Posterior [1746-836]
 Sideless Posterior [1886-211]