



NUMBER OF EVALUATORS
14
COMBINED YEARS IN PRACTICE
272

See what your colleagues have to say about hundreds of other evaluated products

View More Evaluations Online
www.dps.li/p/product-reviews

BENCO DENTAL

e-on Flowable

FLOWABLE COMPOSITE

Benco Dental's e-on Flowable combines nanotechnology and resin-based chemistry to offer optimal flow with the strength of a universal.

An evaluator with 22 years of experience said, "I really liked the product. It is easy to use and flows easily. I have used a lot of flowables over the years and this is one of the best." This evaluator rated all features of e-on Flowable as excellent.

Benco Dental's e-on Flowable is a nanohybrid, light-cured composite that combines innovative nanotechnology and resin-based chemistry. It features an 80% fill, excellent polishability, radiopacity, and low shrinkage.

According to Benco Dental, e-on Flowable offers optimal flow with the strength of a universal composite, and it is suited for both anterior and posterior restorations. In addition, it can be used with all conventional bonding agents and is available in syringe delivery in 7 shades for optimal esthetics.

Fourteen *Dental Product Shopper* evaluators, with an average experience level of more than 14 years in dentistry, used e-on Flowable in their practices for 4 weeks. The evaluators rated this flowable composite on its handling, esthetics, and a number of other features.

Ease of Use

Manufacturer Description: Benco Dental describes e-on Flowable as having excellent handling in a convenient syringe delivery system. It can be used with any standard bonding age, so the clinician can continue to use their preferred adhesive system.



Evaluator Feedback: An evaluator with 34 years of experience said, "The flow and adaptability made for excellent handling." Another evaluator liked the handling and said e-on Flowable was "not as runny as other flowable composites."

An evaluator that practices in Tennessee stated, "The syringe would continue to extrude material that if you did not use it right...requiring you to get another tip." This same evaluator went on to rate his overall satisfaction with e-on Flowable as very good.

An evaluator with 38 years of experience liked best the "excellent consistency and flow" of e-on Flowable, saying it "self-leveled nicely." In addition, he rated all features of the composite as very good or excellent.

Another evaluator said he "enjoyed the honey-like consistency" and added that he would probably purchase e-on Flowable in the future.



The flow and adaptability made for excellent handling.

Eric Van Zytveld, DDS
Denver, CO

4 PREV
EVALS **34** YEARS
EXP



Shade blending could be improved.

Kane Fulton, DMD
Birmingham, AL

16 PREV
EVALS **10** YEARS
EXP




It self-leveled nicely.

William Lieberman, DDS
Red Bank, NJ


12 PREV
EVALS **38** YEARS
EXP

Versatility and Esthetics

Manufacturer Description: Benco Dental states that e-on Flowable offers the strength of a universal composite with low shrinkage and excellent esthetics, making it suitable for both anterior and posterior indications. In addition, it features high radiopacity and easy polishability for convenience.



Evaluator Feedback: An evaluator that practice in Denver, CO, said, “Radiopacity was excellent.” This evaluator called e-on Flowable somewhat better than similar products offered by other manufacturers.



An evaluator with 5 years of experience said, “It seemed strong enough for posterior teeth, but I don’t know if I’d feel comfortable enough placing it routinely in large post fillings.” He went on to rate e-on Flowable as very good overall and said he would probably recommend it to colleagues. Another evaluator called it “a nice product that works well in Class 5 preps.”

One of the evaluators said he would like to see the esthetics improved, saying, “Color stability seemed a little off (too grey).”

Overall Satisfaction

At the end of e-on Flowable evaluation, 11 of the 14 evaluators said they would definitely or probably purchase e-on Flowable in the future, and the same number said they would recommend it to colleagues.

An evaluator who practices in Ohio said, “Initial placement was great with

no observed negative points [compared to my] current flowable.”

When asked what could be improved about e-on Flowable, an evaluator with 18 years of experience said, “All is satisfactory.”

BENCO DENTAL

295 Centerpoint Blvd
Pittston, PA 18640
800.GO.BENCO
www.benco.com

EVALUATION SNAPSHOT

